

Snooooooze
a flexible sleep system

Advent
MEDICAL

Andy Haines

Tel: 07507 850934

E: sales@adventmedicalequipment.co.uk

EARLY INTERVENTION

Suitable and used



in Palliative Care



SNOOOOOOZE

The range is built on the skills of carers and health professionals who have for so long used rolled up towels, pillows, cushions and of course Teddy Bears to position children and we have redesigned them to give cushioning which performs. Either by interfacing to a MATTRESS COVER or for day time, to a THERAPY MAT, the system is easy to use and offers a high degree of user compliance; and changes of posture are very easy with no fixed shapes, offering a high degree of carer compliance. Covered in Dartex, a breathable healthcare material, and filled with crumbed visco-elastic foam there is an excellent tactile feel. All elements are rated at 'Crib 5' fire resistance and the removable covers are washable at 95° for infection control. Our Early Intervention range is a selection of items that allow the carer/therapist to establish and maintain a good postural position for the client and for parents to learn good practice and understand the value of a good sleeping posture.

HORSESHOE

With the little ones low body mass it is important to have cushioning that is adjustable in volume, all snooooooze cushions allow you to add and remove filling. The HORSESHOE is primarily a side layer but has gained a reputation for upper body support either with or without extra support. In supine we can also add a WRAP for greater positioning support.

THE STAR

Having the ability to profile the points the STAR offers a head support that will give the required resting position when the child falls asleep, the points offering some control if needed. Can be used side and supine. With two points under the shoulders we encourage a better head posture.

HEAD PAD

A simpler head support that nestles the head to achieve good mid-line posture, can be used side and supine



The HORSESHOE and STAR, showing alignment of the upper leg and the body in balance. A roll Supporting the back and the front offering something to 'cuddle'.



The HORSESHOE here as a supine positioner, shown on the right with types of supplementary positioning. CURVED WEDGE on the right and TRI-FOLD on the left. The third photo show the WRAP option. The wrap is fixed to the HORSESHOE to achieve even draw on both sides.

In supine or side lying the STAR is a very flexible piece, able to be profiled by returning the points



The HEAD PAD is a simple solution for smaller heads. Like any snooooooze element it can be custom sized.

Snooooooze
a flexible sleep system

POSITIONING • COMFORT • SLEEP



SIDNEY, in the park, on the beach, on holiday, temporarily overnight, now available with a WRAP.



W LEG TROUGH, from a simple reminder of where the legs should lie, to greater support with WRAPS and WEDGES



Above with the wrap in place and supplementary support brackets.



Full body positioning from the BODY ROLL, would usually be supplemented with extra support.



The addition of a wedge to the base can alter the rake.

This element allows the feet to be positioned in bed, often used in conjunction with a W LEG TROUGH



POSITIONING • COMFORT • SLEEP

SIDNEY

The concept was designed by a MUM and developed at Peacocks into what we have today, SIDNEY is designed for multi-use during the day introducing a new posture and the associated benefits, stretching for toys, change of position. Our original MUM wanted to be able to hold her son safely, when on holiday, while she put her make up on!

W LEG TROUGH

A long term favourite in all sizes this element allows the legs to lie in correct alignment and with the WRAP pulled over more security of posture is achieved. Like all elements they can be further supported if needed.

BODY ROLL

Used in side posture in adults this roll offers a similar approach to the HORSESHOE for children where a thinner roll is preferable. Will almost always have supporting units and we recommend CURVED WEDGES.

EXTRA SUPPORT

With a variety of supplementary supports we can add greater support as and only when needed. Either, on other elements or directly onto the body. Right of the BODY ROLL we show the three supplementary types of support. WEDGE AND PILLOW top left is a very ridged unit, top right is the TRI-FOLD where we need similar support but space say for an arm. Below is a CURVED WEDGE which by profile nestles very closely to the body shape, for contact to the body the padded option is preferred with a skin of visco-elastic foam.

FOOT POSITIONING

This new piece offers alignment of the feet, can be custom sized for specific postures, still evolving a set of wedges will be available to change the rake.

MATTRESS COVERS

A fitted sheet with draw string tightening, available in all cot sizes gives us a receptive surface for fixing the elements down.

THERAPY MAT

Used in daytime therapy we have them in sensory rooms, used camping, on trampolines in rebound therapy and very creatively in a sand pit! A new use is to make the system more transportable. If the child stays in various homes the system can easily go with them.